

**I have a Medtronic Defibrillator which is under advisory. Rapid battery depletion may occur.**

**See test instructions on other side.**

**Doctor:** If abnormal test parameters:

**MONITOR IN RESUSC ROOM and call:**

**705.749.8419 (0800-1600 Mon-Fri)**

**705.740.7177 (after hrs)**

**Test Procedure:**

**Apply magnet over device for maximum 5 seconds.**

**If NO AUDIBLE TONE, device has FAILED.**

**Continue to monitor and contact clinic.**

**Peterborough Regional Device Clinic**